

030

Page 2 of 2

Question 1.0  $\rightarrow$  A only  
 Question 2.0  $\rightarrow$  B only  
 Question 3.0  $\rightarrow$  A, B no overlap  
 Question 4.0  $\rightarrow$  no A, B

Uspation  
 megas A833351L  
 A833351O  
 slow 100 pps / full  
 50-2 each mixed  
 ↓  
 25  
 40 slow cool

~~Vagueron~~ - 7<sup>m</sup> along (low point)  
- 2nd + 28% over  
  
- 7<sup>m</sup> HAD  
- 1.5m Wx water  
- 1.5 m down AFB.  
~~Liquid~~  
150  
↓ 96° Lk  
↓ 65°C 15'

James W. Fox  
The City Bank



## GIBCO BRL Custom Primers

### Certificate of Analysis

**Primer 1:**

Primer Name: UBI HSP VER. 1A

Researcher:

Primer Number: A8333C10 (C10)

Primer Length: 69

Sequence (5' to 3'): PAG ACG GCA CGG CAT CTC TGT CGC TGC CTC CAC CGT TGG ACT TGC TCC GCT  
GTC GCC ATC CAG AAA TMolecular Weight  $\mu\text{g}/\mu\text{mole}$ : 21299.2

Millimolar Extinction Coefficient: 678.6

Purity: Desalted

Tm (1 M Na<sup>+</sup>): 96Tm (50 mM Na<sup>+</sup>): 76

% GC: 60

Notes:

 $\mu\text{g per OD}$ : 31.3

nmol per OD: 1.4

OD's: 39.3

 $\mu\text{g's}$ : 1234

nmol: 67

Coupling Eff.: 99%

~5%  
~1%**Primer 2:**

Primer Name: UBI HSP VER.1B

Researcher:

Primer Number: A8333C11 (C11)

Primer Length: 67

Sequence (5' to 3'): PTT TCT GGA TGC CGA CAG CCG AGC AAG TCC AAC GGT GGA GGC AGC GAC AGA  
GAT GCC GTG CCG TCT GCMolecular Weight  $\mu\text{g}/\mu\text{mole}$ : 21597.4

Millimolar Extinction Coefficient: 732.6

Purity: Desalted

Tm (1 M Na<sup>+</sup>): 97Tm (50 mM Na<sup>+</sup>): 78

% GC: 62

Notes:

 $\mu\text{g per OD}$ : 29.8

nmol per OD: 1.3

OD's: 10.7

 $\mu\text{g's}$ : 319

nmol: 14

Coupling Eff.: 98%

VER 1A

57 nmol

570  $\mu\text{L}$   $\rightarrow$  100  $\mu\text{mol}/\mu\text{L}$ 

14 nmol

140  $\mu\text{L}$   $\rightarrow$  100  $\mu\text{mol}/\mu\text{L}$ 

\* See Note about Quantities in  
Supporting Information.

LIFE  TECHNOLOGIES.

# GIBCO BRL Custom Primers Certificate of Analysis

**Primer 1:**

Primer Name: UBI HSPA VER.2A

Researcher:

Primer Number: D0373887 (B07)

Primer Length: 81

Sequence (5' to 3'): P-A GAC GGC ACG GCA TCT CTG TCG CTG CCT CTG GAC CGG TGT CGA CCA CCG  
TTG GAC TTG CTC CGC TGT CGG CAT CCA GAA ATMolecular Weight  $\mu\text{g}/\mu\text{mole}$ : 26105.2 $\mu\text{g}$  per OD: 31.5

Millimolar Extinction Coefficient: 824.3

nmole per OD: 1.2

Purity: Desalt

OD's: 90.0

Tm (1 M Na+): 98

 $\mu\text{g's}^*$ : 2850

Tm (50 mM Na+): 77

nmole: 108

% GC: 61

Coupling Eff.: 98%

Notes:

100  $\mu\text{g}$  in 100  $\mu\text{L}$  of water**Primer 2:**

Primer Name: UBI HSPB VER.2B

Researcher:

Primer Number: D0373888 (B08)

Primer Length: 82

Sequence (5' to 3'): P-T TTC TGG ATG CCG ACA GCG GAG CAA GTC CAA CGG TGG TCG AGA GGG GTC  
CAG AGG CAG CGA CAG AGA TGG CGT GCG GTC TGCMolecular Weight  $\mu\text{g}/\mu\text{mole}$ : 26872.4 $\mu\text{g}$  per OD: 29.7

Millimolar Extinction Coefficient: 902.2

nmole per OD: 1.1

Purity: Desalt

OD's: 77.0

Tm (1 M Na+): 99

 $\mu\text{g's}^*$ : 2294

Tm (50 mM Na+): 78

nmole: 86

% GC: 63

Coupling Eff.: 98%

Notes:

100  $\mu\text{g}$  in 100  $\mu\text{L}$  of water**Primer 3:**

Primer Name: UBI HSPA VER.3A

Researcher:

Primer Number: D0372809 (B09)

Primer Length: 81

Sequence (5' to 3'): P-A GAC GGC ACG GCA TCT CTG TCG CTG CCT CTC GAG AGT TCC GCT CCA CCG  
TTG GAC TTG CTC CGC TGT CGG CAT CCA GAA ATMolecular Weight  $\mu\text{g}/\mu\text{mole}$ : 26160.2 $\mu\text{g}$  per OD: 31.5

Millimolar Extinction Coefficient: 830.6

nmole per OD: 1.2

Purity: Desalt

OD's: 88.7

Tm (1 M Na+): 98

 $\mu\text{g's}^*$ : 2783

Tm (50 mM Na+): 76

nmole: 106

% GC: 60

Coupling Eff.: 98%

Notes:

100  $\mu\text{g}$  in 100  $\mu\text{L}$  of water\* See Note about Quantities in  
Supporting Information.

# GIBCO BRL Custom Primers Certificate of Analysis

**Primer 4:**

Primer Name: UBI HSPB VER.3B

Researcher:

Sequence (5' to 3'): P-T TTC TGG ATG CCG ACA GCG GAG CAA GTC CAA CGG TGG AGC GGA ACT CTC  
 GAG AGG CAG CGA CAG AGA TGC CGT GCC GTC TGC

Primer Number: D0373B10 (B10)

Primer Length: 82

Molecular Weight  $\mu\text{g}/\mu\text{mole}$ : 26518.4

Millimolar Extinction Coefficient: 901.3

 $\mu\text{g}$  per OD: 29.7

nmoles per OD: 1.1

Purity: **Desalt**Tm (1 M Na<sup>+</sup>): 99Tm (50 mM Na<sup>+</sup>): 77

% GC: 62

Notes:

OD's: 81.2

 $\mu\text{g/s}$ : 2478

nmoles: 32

Coupling Eff.: 98%

930  $\mu\text{g}$   $\rightarrow$  100  $\mu\text{g}$ **Primer 5:**

Primer Name: UBI HSPA VER.4A

Researcher:

Sequence (5' to 3'): P-A GAC GGC ACG GCA TCT CTG TCG CTG CCT CTG GAC CCG TCT CGA CTC GAG  
 AGT TCC GCT CCA CCG TTG GAC TTG CTC CGC TGT CCG CAT CCA GAA AT

Primer Number: D0373B11 (B11)

Primer Length: 86

Molecular Weight  $\mu\text{g}/\mu\text{mole}$ : 30988.2

Millimolar Extinction Coefficient: 978.3

 $\mu\text{g}$  per OD: 31.7

nmoles per OD: 1.0

Purity: **Desalt**Tm (1 M Na<sup>+</sup>): 100Tm (50 mM Na<sup>+</sup>): 78

% GC: 61

Notes:

OD's: 89.3

 $\mu\text{g/s}$ : 2833

nmoles: 91

Coupling Eff.: 98%

930  $\mu\text{g}$   $\rightarrow$  100  $\mu\text{g}$ **Primer 6:**

Primer Name: UBI HSPB VER.4B

Researcher:

Sequence (5' to 3'): P-T TTT TGG ATG CCG ACA GCG GAG CAA GTC CAA CGG TGG AGC GGA ACT CTC  
 GAG TCG AGA CCG GTC CAG AGG CAG CGA CAG AGA TGC CGT GCC GTC TGC

Primer Number: D0373B12 (B12)

Primer Length: 87

Molecular Weight  $\mu\text{g}/\mu\text{mole}$ : 31781.4

Millimolar Extinction Coefficient: 1070.6

 $\mu\text{g}$  per OD: 29.6

nmoles per OD: 0.9

Purity: **Desalt**Tm (1 M Na<sup>+</sup>): 100Tm (50 mM Na<sup>+</sup>): 79

% GC: 62

Notes:

OD's: 97.1

 $\mu\text{g/s}$ : 2883

nmoles: 90

Coupling Eff.: 98%

930  $\mu\text{g}$   $\rightarrow$  100  $\mu\text{g}$ 

\* See Note about Quantities in  
 Supporting Information.



LIFE  TECHNOLOGIES.

9-007

Do Nucleobond prep of 5596, 5597,  
4216, 4217, 4218 and 4219

Digest 5596 and 5597 w/ EcoRI as a check  
make sure smaller frag. is ~ 2kb.

Digest 4216, 4217, 4218, 4219 w/ BglI/SalI  
to check that the 168bp frag. is generated.

Digest 4216 w/ BglI/XbaI to use as  
accepting vector for ubiquitin versions 1-4.  
Gel isolate on 10% agarose

Digest 5596 w/ NheI/NotI to isolate insert  
(~ 1.6 kb) L4-BASS. NotI and NheI  
Gel isolate on agarose. cut in different buffer  
cut w/ NotI 1st. Debris  
and EtOH ppt.  
cut w/ NheI

Mini-preps on BASS:NA #5 5648 - 5665  
Digest w/ EcoRI/PstI. Cut 118 also.  
Run on 10% agarose.


→ Anneal Ubiquitin versions 1, 2, 3, and 4  
aligns together.  
Heat to 95°C for 5 min then stick on ice  
check these on a 10% agarose gel.

Run pre-cut 3270 NheI/NotI on agarose  
gel to check it out. ~~rather simple~~

13-004

Check the variation w/ a 113/110T digest (first 5 of each)

1. 1 Kb ladder  
 2-6. Lbri 1 7543 - 7597  
 7-11. Lbri 2 7561 - 7565  
 12-16. Lbri 3 7579 - 7583  
 17-21. Lbri 4 7597 - 7601  
 22. ~~5000~~ NL/110T 7062 vector  
 23. 1939 NL/110T Spg  
 24. 1 Kb ladder 0




App look fine. Send  
 these for sequencing

Sequencing data shows 7368 to be correct for  
 G43/P2P

Make Avidin/P2P using 7368 as standard

Digest 7424 (Avidin/3220) and 7368 (G43/P2P)  
 with BamHI/HaeIII

Gels: 1. 1 Kb ladder  
 2-5. 7368 (G43/P2P) Bam/Hae  
 6. Spg  
 7-10. 7424 (Avidin) Bam/Hae



Isolate large vector frag from G43/P2P and the smaller  
 insert band from Avidin

Katherine Harper



Gus Assay

PURPOSE: TO QUANTITATE THE AMOUNT OF GUS IN CORN SEED

MATERIALS: REACTION PLATE - COSTAR EIA/RIA  
 READING PLATE - NUNC FLUORENCE POLYSOEP  
 MU - 4 METHYLUMLBELLIFERONE (SIGMA H-1508)  
 MUG - 4 METHYLUMLBELLIFERONE B GLUCURONIDE (SIGMA H-9130)  
 MICROPLATE  
 FLUORESCENCE MICROPLATE READER

PROCEDURE: USE PROTOCOL FOUND ON PAGE # 57 OF THIS NOTEBOOK (#58)

RESULTS: DATA FOUND BELOW (BASED ON 20-NEW READINGS)

SAMPLE #	%TSP	SAMPLE #	%TSP
05E12020-4	0.088	05E 05230-1	0.087
-5	ND	-2	0.54
05D 01120-1	ND	-3	0.61
-2	ND	-4	0.16
-3	ND	-5	0.06
-4	ND	11 0828-1	0.021
-5	ND	-2	0.002
05E 15070-4	0.028	-3	0.007
11 05050-1	0.17	-4	ND
-2	0.015	-5	0.001
-3	0.010	11 07050-1	0.3
-4	0.174	-2	0.089
-5	0.010	-3	0.27
11 05090-1	0.043	-4	0.013
-2	0.014	-5	0.43
-3	0.001		
-4	0.004		
-5	0.004		
05D 01010-1	0.006		
-2	0.010		
-3	0.009		
-4	0.60		
-5	0.048		

Investigator:

Book # 58

Chris Brook Date:

Witness:

Elizabeth Wilcox Date:



GUS ASSAY

SEE PURPOSE, MATERIALS, &amp; PROCEDURE BELOW.

Purpose: To quantify the amount of GUS in corn seed extracts.

Materials: Reaction Plate-96 Wells, 96-well bottom, 96-well flat bottom plate  
 Reaction Plate-96 Wells, 96-well black plate  
 MUG 4-methylumbelliferyl-β-D-glucuronide (Sigma M-1500)  
 MUG 4-methylumbelliferyl-β-D-glucuronide (Sigma M-1510)  
 Microplate reader (Spectrophotometer)

Reagents: 1.0 M Tris-HCl, 50 mM sodium phosphate, pH 7.0, 1 mM EDTA, 10 mM MgCl<sub>2</sub>

Note: 50 mM sodium phosphate is made by mixing 99 ml of stock A (0.2 M Na<sub>2</sub>HPO<sub>4</sub>, 0.2 M NaH<sub>2</sub>PO<sub>4</sub>) with 1 ml of stock B (0.2 M NaOH, 0.2 M NaCl) and bringing to a final volume of 1.0 L with H<sub>2</sub>O.

Also note that the 10 mM EDTA should be added to an aliquot of the Tris buffer from daily, enough for that day's experiment.

Stock Buffers: 0.2 M Na<sub>2</sub>HPO<sub>4</sub> (0.2 M)  
 1 M NaCl (1 M)  
 0.2 M NaOH (0.2 M)  
 0.2 M NaCl (0.2 M)

Procedure: Corn seed extracts should always be prepared according to the total protein assay using the standard procedure.

In a reaction plate, add 100 μl of seed extract to a total volume of 100 μl Tris buffer. GUS activity can be assayed with 1 μg total protein. Samples should be assayed in triplicate.

Add standard curve to triplicate wells as follows:

10 μl of 1.0 M Tris-HCl standard stock to diluted with 90 μl Tris buffer, 10 μl of this 1:10 dilution to further diluted with 90 μl Tris buffer to give a 1:100 dilution.

0.1 M Tris-HCl standard  
 1000 μl Tris-HCl standard  
 10,000 μl Tris-HCl standard  
 100,000 μl Tris-HCl standard

100 μl Tris buffer / well  
 12.5 μl of the 1:100 dilution = 12.5 μl Tris buffer / well  
 12.5 μl of the 1:10 dilution = 12.5 μl Tris buffer / well  
 12.5 μl of the 1:1000 dilution = 12.5 μl Tris buffer / well

CB

Prepare the reaction plates by pipetting 175 μl of Tris buffer into each well of the plate. You will need a separate plate for each time point required. Generally we take readings at 0, 15, 30 and 45 minutes.

Then the 20 μl MUG substrate stock to 3 wells with Tris buffer. Add 25 μl of 5 mM MUG to every well including both standard and sample wells and mix to start the reaction. Immediately after adding the MUG, pipette 25 μl of seed extract from the reaction plate into a prepared reaction plate. Then the reaction plate is at 37 °C until the next time point. At each subsequent time point, pipette 25 μl of seed extract from the reaction plate into a prepared reaction plate.

Reaction is stable for several hours once it has been stopped. Note that stopping the reaction is essential for fluorescence detection.

Plates are read at 540 nm excitation wavelength and 490 nm emission wavelength.

The unknown samples are read against the standard curve to get MUG and the amount of GUS in the sample is calculated as follows:

Assume that MUG is the only substrate (from Value Column) of substrate reaction provided = 0.1 M MUG / ml \* 60 min / hr = 0.1 MUG / hr. Note that if there is a 1:100 dilution of the sample, then the value must be multiplied from the average MUG of each subsequent reading. This value is then converted to the amount of protein added to the sample by dividing by the total protein added to give μg MUG / hr / μg. This value is converted to GUS by multiplying by 1.0E+10 which is a conversion factor determined value of Fluorescence.

A Quality Control sample (a known amount of GUS) is added to each assay to determine reproducibility of quantitation.

CB

## RESULTS: DATA FOUND BELOW. (10-MIN READINGS)

Sample#	%TSP	Sample#	%TSP	Sample#	%TSP
GSC 01040-1	0.0	GSC 01110-1	0.0	GSC 01060-1	0.0
-2	0.4 0.04	-2	0.4 0.04	-2	0.0
-3	0.6 0.06	-3	0.0	-3	0.0
-4	0.5 0.05	-4	0.4 0.04	-4	0.0
-5	0.4 0.04	-5	0.4 0.04	-5	4.8 0.5
GSD 02130-1	1.1 0.1	GSC 01070-1	4.2 0.4	GSC 01130-1	8.4 0.8
-2	0.7 0.07	-2	2.7 0.3	-2	0.1 0.01
-3	0.9 0.1	-3	3.7 0.3	-3	8.6 0.9
-4	0.0	-4	5.2 0.5	-4	5.0 0.5
-5	0.8 0.1	-5	0.01 0.001	-5	0.7 0.07
GSC 01020-1	0.0	GSC 01040-1	0.1 0.01	GSC 01110-1	0.0
-2	0.0	-2	5.1 0.5	-2	9.2 0.9
-3	0.12 0.01	-3	0.3 0.03	-3	0.0
-4	0.0	-4	0.3 0.03	-4	0.0
-5	0.2 0.02	-5	0.04 0.004	-5	9.6 0.7
GSC 01030-1	0.0				
-2	4.0 0.4				
-3	4.2 0.4				
-4	0.5 0.05				
-5	9.5 0.8				

Investigator:

Investigator: Chris Bivola

Block # 167

Date:

Witness:

Witness: Elizabeth Wilfong

Date: